


SONIC



1. Hold the respirator in your hand with the nose piece at your fingertips. Allow the headbands to hang freely below your hand.

2. Cup the respirator firmly against your face with the nose piece on the bridge of your nose.

3. Stretch and pull the lower headband and position below your ears. Stretch and pull the top headband on the back of your head and above your ears.

4. Press and form the soft nose piece to the shape of your nose.

5. To test, cup both hands over the respirator and exhale vigorously.

- If airflows around your nose, tighten the nose piece.
- If air leaks around the edges, reposition the headband for a better fit.